BECKMAN COULTER

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Optima L / XL / LE serial Operation

STANDARD OPERATION

1. With the power on, open the chamber door, install the rotor, and close the chamber door. Press [VACUUM] if you want to evacuate the chamber now.
2. Enter the required run conditions (SPEED, TEMP, and TIME, HOLD, or $\omega 2t$).
3. Select ACCEL and DECEL profiles.
4. Press [ENTER/RECALL] and [START]. (The vacuum system is automatically activated now, unless you pressed [VACUUM] earlier.)

Press [STOP] to terminate a run in the HOLD mode (or to stop any run for any reason). Runs in the timed and $\omega 2t$ modes will terminate automatically and the instrument will sound a tone to indicate the end of the run. After the rotor has stopped, press [VACUUM] to vent the chamber.

Open the chamber door and remove the rotor. To keep the rotor chamber clean and dry, keep the door closed between runs.

REPEATING A RUN

To repeat a run, press [ENTER/RECALL] and [START]. There is no need to reenter run conditions unless you wish to make a change.

PROGRAMMED OPERATION

1. With the power on, open the chamber door, install the rotor, and close the chamber door. Press [VACUUM] if you want to evacuate the chamber now.
 2. Press [PROG] and use the keypad to select a program number.
(Or press [PROG] repeatedly until the required program appears.)
- 3-3 Run Procedure
3. Press [ENTER/RECALL] and [START] to begin the run. (The vacuum system is automatically activated now, unless you pressed [VACUUM] earlier.)

Any of the run conditions can be changed while a programmed run is in progress. Such changes will affect only the current run. To change a parameter during the

run:

1. Press a parameter key (for example, [SPEED]).
2. Enter the new value.
3. Press [ENTER/RECALL].

The run in progress will be modified accordingly. Runs in the timed and ω 2t modes will terminate automatically when the set value is reached.

4. After the rotor has stopped, press [VACUUM] to vent the chamber.
5. Open the chamber door, remove the rotor, and close the door.

To keep the rotor chamber clean and dry, keep the door closed between runs.

Maintenance and Clean:

1. Keep Chamber Clean
2. Need check spinning speed reached user's setup speed rpm
3. Keep Rotor Clean and Rotor Lid O-ring need clean with vacuum grease °
4. Keep Door O-ring Clean with vacuum grease °